


BEDFORD
ElastiTab®
 Simple on-pack sampling tool




ElastiTab® is a versatile option for carrying sachets, packets, or other lightweight items. The combination of a tab with pressure sensitive adhesive, and a stretchy elastomer loop—gives you an easy-to-apply device for on-pack promotions!




ElastiTab® standard options:




Die #2102
Recommended for items with a circumference of 4.4"–6.5" (111mm–165mm)




Die #2144
Recommended for items with a circumference of 2.7"–4.5" (68mm–114mm)




Die #2145
Recommended for items with a circumference of 4.5"–8.2" (114mm–208mm)




Die #2243
Recommended for items with a circumference of 1.59"–3.33" (40mm–84mm)



Die #2164
Recommended for items with a circumference of 4.5"–8.2" (114mm–208mm)



Die #2347
Recommended for items with a circumference of 2.7"–4.5" (68mm–114mm)



Die #2312
Recommended for items with a circumference of 2.6"–3.8" (66mm–97mm)

Details:

- ▶ A 2 inch ElastiTab® with pressure sensitive adhesive has the ability to hold up to a 2 oz. packet/sachet.
- ▶ A 1.5 inch or 1 inch ElastiTab® with pressure sensitive adhesive has the ability to hold up to a 1 oz. packet/sachet.
- ▶ For ease of application, the ElastiTab® can be hand or machine applied to your packet/sachet. Call us for machine information.
- ▶ Elastomer loops are available in a variety of colors.

If you're looking for a custom option, we would love to work on an alternative solution together.

We will need to know...

- ▶ The weight, size, and material of the item you would like to carry with the ElastiTab®.
- ▶ The circumference of the item where the elastomer loop will be applied.



Bedford Industries | ElastiTag®

1659 Rowe Ave | PO Box 39
Worthington, MN 56187-0039
Ph: 1-507-376-4136 | Fax: 1-507-376-6742
Toll Free US/Canada:
1-877-BEDFORD (233-3673)
www.bedford.com | www.elastitag.com
email: bedford@bedford.com

Contact Us Today!